


LASCHAL[®] SURGICAL

ENDOTECH NZ CATALOGUE

Laschal endodontic instruments for control, access, and retrieval.

A clinician-facing catalogue for tactile file control, separated instrument retrieval, gutta percha removal, and microsurgical endodontic support.



83

PRODUCTS
LISTED

5

FEATURED
PATHWAYS

13

VIDEO
LINKS

CLINICAL WORKFLOW FIT

Choose by the clinical problem in front of you.

Laschal instruments are most useful where direct finger access, visibility, or force control limits the procedure.

01

File control and pre-curving

Bend, place, and manipulate files with less finger crowding when anatomy or microscope angle limits direct access.

02

Separated instrument retrieval

Use angled probes and Steiglitz forceps as a sequence for exposing, troughing, and removing fragments.

03

Retreatment and gutta percha removal

Keep retreatment instruments organized by direction and movement so removal remains visible and controlled.


04

Microsurgical access and closure

Needle holders, scissors, tissue forceps, and suture support instruments round out the endodontic surgical setup.

FEATURED INSTRUMENTS

Core Laschal pathways for endodontic access and retrieval.



PART NUMBER: 90AHF/L, 75CHF/L, EF-1-7

Tactile endo file forceps

CLINICAL PROBLEM

Deep posterior access and microscope positioning can make small files difficult to control with fingers alone.

LASCHAL SOLUTION

Angled file forceps place file handling into a slim instrument profile for tactile access.

KEY FEATURES

- 75 or 90 degree access options
- Diamond-coated working surfaces
- Tactile file placement support

CLINICAL BENEFITS

- Improves file control under magnification
- Reduces finger crowding around the access
- Supports controlled placement in difficult anatomy

TECHNIQUE VIDEOS

[Watch Dr. Marc Habib use EF-1-7 tactile file forceps](https://www.youtube.com/watch?v=O14DwrteW9M)

<https://www.youtube.com/watch?v=O14DwrteW9M>

[Watch Laschal endodontic file forceps overview](https://www.youtube.com/watch?v=D0TvF82CMR4)

<https://www.youtube.com/watch?v=D0TvF82CMR4>



PART NUMBER: LA-F

Endodontic file bender

CLINICAL PROBLEM

Negotiating complex canal anatomy often requires a repeatable pre-curve rather than a rough hand bend.

LASCHAL SOLUTION

LA-F is the current Laschal file-bending tool for controlled endodontic file pre-curving.

KEY FEATURES


- Straight Uniband control
- File bending and pre-curving support
- Compact chairside format

CLINICAL BENEFITS

- Helps preserve intended file path
- Supports controlled negotiation
- Pairs naturally with tactile file forceps

FEATURED INSTRUMENTS

Core Laschal pathways continued.



PART NUMBER: FXP/S, FXP/45, FXP/75, FXP/90, FXP/110

FXP file extraction protocol

CLINICAL PROBLEM

Separated instruments can require a narrow, controlled retrieval pathway without over-enlarging dentine.

LASCHAL SOLUTION

Angled diamond excavation probes and Steiglitz forceps create a staged pathway for ultrasonic troughing and retrieval.

KEY FEATURES

- 45, 75, 90, and 110 degree probe options
- Probe bender support
- Steiglitz forceps pairing

CLINICAL BENEFITS

- Supports a conservative retrieval strategy
- Keeps access choices visible
- Gives clinicians a repeatable protocol sequence

TECHNIQUE VIDEOS

Watch the FXP troughing system


<https://www.youtube.com/watch?v=afvIQpxp7Ls>

Watch Dr Ortiz Hughes demonstrate FXP retrieval

<https://www.youtube.com/watch?v=249Xa1SXEd8>

Watch Steiglitz forceps retrieval support

<https://www.youtube.com/watch?v=iJ9l7UQPtVs>



PART NUMBER: SET-GP

Gutta percha removal set

CLINICAL PROBLEM

Retreatment requires removal of gutta percha while preserving visibility and avoiding unnecessary canal disruption.

LASCHAL SOLUTION

A directional removal set uses scraper-style instruments that can be pre-heated for controlled gutta percha removal.

KEY FEATURES

- Four directional instruments
- Sterile tray organization
- Optional heat support

CLINICAL BENEFITS

- Keeps the working field organized
- Supports controlled scraping movement
- Fits a retreatment-focused EndoTech workflow

FEATURED INSTRUMENTS

Core Laschal pathways continued.



PART NUMBER: 45-S/L, 75-SP/L

Raptor post and root-tip forceps

CLINICAL PROBLEM

Posts, points, and small fragments often need narrow access with force delivered only where it is useful.

LASCHAL SOLUTION

Raptor forceps use a slim, durable tip profile for controlled gripping in tight endodontic and surgical spaces.

KEY FEATURES

- 45 and 75 degree options
- N/S and E/W orientation choices
- Narrow carbide tip geometry

CLINICAL BENEFITS

- Improves approach options
- Supports low-force grip control
- Extends the endodontic set beyond canal instrumentation

FULL ENDODONTIC RANGE

File retrieval

1-6 OF 18 PRODUCTS



FILE RETRIEVAL

15cm Universal RAPTOR Forceps for points, posts and pedodontic crowns ThumLok® Carbide Flat Handle 45 deg N/S, narrow and durable Tips

SKU: 45-S/L

Enquire through EndoTech NZ



FILE RETRIEVAL

15cm Micro Steiglitz Forceps ThumLok® Carbide Flat Handle 75 deg E/W

SKU: 75SPL/M

Enquire through EndoTech NZ



FILE RETRIEVAL

14.75cm Universal RAPTOR Forceps ThumLok® Carbide Flat Handle 75 deg E/W

SKU: 75-SP/L

Enquire through EndoTech NZ



FILE RETRIEVAL

14.75cm Universal RAPTOR Forceps Carbide Flat Handle 75 deg E/W

SKU: 75-SP

Enquire through EndoTech NZ



FILE RETRIEVAL

14.5cm Micro Steiglitz Forceps ThumLok® Carbide Flat Handle 75 deg N/S

SKU: 75SL/M

Enquire through EndoTech NZ



FILE RETRIEVAL

14.25cm Universal Steiglitz Forceps ThumLok® Carbide Flat Handle 90 deg N/S

SKU: 90SL/M

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

File retrieval continued

7-12 OF 18 PRODUCTS



FILE RETRIEVAL

14.5cm Universal Steiglitz Forceps ThumLok® Carbide Flat Handle 90 deg E/W

SKU: 90SPL/M

Enquire through EndoTech NZ



FILE RETRIEVAL

15cm Micro Steiglitz Forceps ThumLok® Diamond 45 deg N/S Angulation

SKU: D-45SL/M

Enquire through EndoTech NZ




FILE RETRIEVAL

16cm Micro Steiglitz Forceps ThumLok® Diamond 45 deg E/W Angulation

SKU: D-45SPL/M

Enquire through EndoTech NZ



FILE RETRIEVAL

15.5cm Micro Steiglitz Forceps ThumLok® Diamond 90 deg E/W Angulation

SKU: D-90SPL/M

Enquire through EndoTech NZ




FILE RETRIEVAL

14.25cm Micro Steiglitz Forceps ThumLok® Diamond 90 deg N/S Angulation

SKU: D-90SL/M

Enquire through EndoTech NZ



FILE RETRIEVAL

15cm Micro Steiglitz Forceps ThumLok® Diamond 75 deg E/W Angulation

SKU: D-75SPL/M

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

File retrieval continued

13-18 OF 18 PRODUCTS



FILE RETRIEVAL

**15cm Micro Steiglitz Forceps ThumLok®
Diamond 75 deg N/S Angulation**

SKU: D-75SL/M

Enquire through EndoTech NZ



FILE RETRIEVAL

Laschal File Extraction Protocol

SKU: FXP/S

Enquire through EndoTech NZ




FILE RETRIEVAL

11.6cm File Extraction Probe 90 deg

SKU: FXP/90

Enquire through EndoTech NZ



FILE RETRIEVAL

12.5cm File Extraction Probe 75 deg

SKU: FXP/75

Enquire through EndoTech NZ



FILE RETRIEVAL

11.6cm File Extraction Probe 110 deg

SKU: FXP/110

Enquire through EndoTech NZ



FILE RETRIEVAL

14cm File Extraction Probe 45 deg

SKU: FXP/45

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Endo file control

4 PRODUCTS



ENDO FILE CONTROL

**15cm Universal Endo File Forceps ThumLok®
Diamond Flat Handle 75 deg E/W**

SKU: 75CHF/L

Enquire through EndoTech NZ



ENDO FILE CONTROL

**15.75cm Universal Endo File Forceps
ThumLok® Diamond Flat Handle 90 deg N/S**

SKU: 90AHF/L

Enquire through EndoTech NZ



ENDO FILE CONTROL

**17.75cm Universal Endo File Forceps
ThumLok® Diamond Round Handle Straight
Straight Blunt Tips**

SKU: EF-1-7

Enquire through EndoTech NZ



ENDO FILE CONTROL

**11.5cm Uniband Forceps Straight Standard
Tips**

SKU: LA-F

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Needle holders

1-6 OF 25 PRODUCTS



NEEDLE HOLDERS

14cm Silk 3-0/7-0, CGS 4/0-7/0 Webster Needle Holder Ring-Handled Carbide Straight

SKU: 2-331-12

Enquire through EndoTech NZ



NEEDLE HOLDERS

16.5cm Silk 2-0/6-0, CGS 3-0/7-0 Needle Holder Ring-Handled Carbide Straight

SKU: 2-331-15

Enquire through EndoTech NZ



NEEDLE HOLDERS

16.5cm Silk 1/4-0, CGS 0/5-0 Mayo-Hegar Needle Holder with Side Cutting Suture Scissors Ring-Handled Carbide Straight

SKU: 2-331-20

Enquire through EndoTech NZ



NEEDLE HOLDERS

14cm Silk 3-0/7-0, CGS 4/0-7/0 Needle Holder Ring-Handled Carbide Straight

SKU: 2-332-12

Enquire through EndoTech NZ




NEEDLE HOLDERS

14cm Silk 1/4-0, CGS 0/5-0 Needle Holder Ring-Handled Carbide Straight

SKU: 2-332-20

Enquire through EndoTech NZ



NEEDLE HOLDERS

17.5cm Silk 1/4-0, CGS 0/5-1 Needle Holder Ring-Handled Carbide Straight

SKU: 2-335-20

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Needle holders continued

7-12 OF 25 PRODUCTS



NEEDLE HOLDERS

17.75cm Universal Cutting-Edge Needle Holder with Suture Cutter Carbide Round Handle Straight Tips

SKU: 7-10RL

Enquire through EndoTech NZ



NEEDLE HOLDERS

15cm Universal Cutting-Edge Needle Holder with Suture Cutter ThumLok® Carbide Round Handle Curved Tips

SKU: 6-10RLC/TL

Enquire through EndoTech NZ



NEEDLE HOLDERS

15cm Universal Cutting-Edge Needle Holder with Suture Cutter ThumLok® Carbide Round Handle Straight Tips

SKU: 6-10RL/TL

Enquire through EndoTech NZ



NEEDLE HOLDERS

15cm Universal Cutting-Edge Needle Holder with Suture Cutter Carbide Round Handle Curved Tips

SKU: 6-10RL/C

Enquire through EndoTech NZ



NEEDLE HOLDERS

15cm Universal Cutting-Edge Needle Holder with Suture Cutter Carbide Round Handle Straight Tips

SKU: 6-10RL

Enquire through EndoTech NZ



NEEDLE HOLDERS

17.75cm Universal Cutting-Edge Needle Holder with Suture Cutter ThumLok® Carbide Round Handle Curved Tips

SKU: 7-10RLC/TL

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Needle holders continued

13-18 OF 25 PRODUCTS



NEEDLE HOLDERS

**17.75cm Universal Needle Holder
Castroviejo Lock Carbide Round Handle
Straight Tips**

SKU: 7-TCL/R

Enquire through EndoTech NZ



NEEDLE HOLDERS

**17.75cm Universal Cutting-Edge Needle
Holder with Suture Cutter Carbide Round
Handle Curved Tips**

SKU: 7-10RL/C

Enquire through EndoTech NZ



NEEDLE HOLDERS

**18cm Universal Needle Holder ThumLok®
Carbide Round Handle WIDE JAWS Curved
Tips**

SKU: 7-TCLCRW/TL

Enquire through EndoTech NZ



NEEDLE HOLDERS

**18cm Universal Needle Holder Castroviejo
Lock Carbide Round Handle WIDE Tips**

SKU: 7-TCLR/W

Enquire through EndoTech NZ



NEEDLE HOLDERS

**18cm Universal Needle Holder ThumLok®
Carbide Round Handle Straight Tips**

SKU: 7-TCLR/TL

Enquire through EndoTech NZ



NEEDLE HOLDERS

**17.75cm Micro Needle Holder ThumLok®
Carbide Round Handle Straight Micro Tips**

SKU: 7TCLR/MM/TL

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Needle holders continued

19-24 OF 25 PRODUCTS



NEEDLE HOLDERS

17.75cm Micro Needle Holder Castroviejo Lock Carbide Round Handle Straight Micro Tips

SKU: 7TCLR/MM

Enquire through EndoTech NZ



NEEDLE HOLDERS

18cm Universal Needle Holder Castroviejo Lock Carbide Flat Handle WIDE JAWS Curved Tips

SKU: 7-TCLR/W

Enquire through EndoTech NZ



NEEDLE HOLDERS

18cm Universal Needle Holder ThumLok® Carbide Round Handle Curved Tips

SKU: 7-TCLR/TL

Enquire through EndoTech NZ



NEEDLE HOLDERS

17.5cm Micro Needle Holder ThumLok® Carbide Round Handle Curved Tips

SKU: 7TCLR/MM/TL

Enquire through EndoTech NZ



NEEDLE HOLDERS

17.5cm Micro Needle Holder Castroviejo Lock Carbide Round Handle Curved Tips

SKU: 7TCLR/MM

Enquire through EndoTech NZ



NEEDLE HOLDERS

18cm Universal Needle Holder Castroviejo Lock Carbide Round Handle Curved Tips

SKU: 7-TCLC/R

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Needle holders continued

25-25 OF 25 PRODUCTS



NEEDLE HOLDERS

**18cm Universal Needle Holder ThumLok®
Carbide Round Handle WIDE JAWS Straight
Tips**

SKU: 7-TCLRW/TL

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Microsurgical scissors

1-6 OF 21 PRODUCTS



MICROSURGICAL SCISSORS

13.5cm Periodontal Scissors Castroviejo Flat Handle 45 deg angle 1.0 cm Curved Blade

SKU: 51-12-45C

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

14.5cm Periodontal Scissors Castroviejo Flat Handle 45 deg angle 1.75 cm Curved Blade

SKU: 51-15-45C

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

14.5cm Periodontal Scissors Castroviejo Flat Handle 30 deg angle 1.75 cm Curved Blade

SKU: 51-15-30C

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

13.5cm Periodontal Scissors Castroviejo Flat Handle 30 deg angle 1.0 cm Curved Blade

SKU: 51-12-30C

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

13.5cm Vannas Scissors Castroviejo Flat-Handle Straight 1cm Straight Sharp Tips

SKU: DS-1

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

14cm Reach Anywhere Scissors Castroviejo Flat-Handle Straight 1.25 cm Straight Sharp Tips

SKU: DS-G

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Microsurgical scissors continued

**7-12 OF 21
PRODUCTS**



MICROSURGICAL SCISSORS

**9cm Vannas Scissors Castroviejo Flat-Handle
Straight 3 mm Straight Sharp Tips**

SKU: DS-1M

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**9cm Vannas Scissors Castroviejo Flat-Handle
3mm Curved Sharp Tips**

SKU: DS-1CM

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**13.5cm Vannas Scissors Castroviejo Flat-
Handle 1cm Curved Sharp Tips**

SKU: DS-1C

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**13cm Scissors/Forceps Combination
Castroviejo Flat Handle Angled Sharp Tips**

SKU: N-103F

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**18cm Scissors Castroviejo Round Handle
Curved Sharp Tips**

SKU: N-1C-7

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**14cm Scissors Castroviejo Flat Handle
Curved Sharp Tips**

SKU: N-1C

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Microsurgical scissors continued

**13-18 OF 21
PRODUCTS**



MICROSURGICAL SCISSORS

**17.5cm Scissors Castroviejo Round Handle
Straight Sharp Tips**

SKU: N-1-7

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**14cm Scissors Castroviejo Flat Handle
Straight Sharp Tips**

SKU: N-1

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**13cm Scissors Castroviejo Flat Handle
Straight Sharp Tips**

SKU: N-4

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**13.5cm Littauer Scissors Castroviejo Flat
Handle 45 deg N/S Straight Sharp Tips**

SKU: N-6A

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**13.5cm Littauer Scissors Castroviejo Flat
Handle Straight Sharp Tips**

SKU: N-6

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**SofTouch Suture Scissors 13.5cm Castroviejo
Flat Handle Straight Blunt Tips Duck Billed**

SKU: N-4XF

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Microsurgical scissors continued

**19-21 OF 21
PRODUCTS**



MICROSURGICAL SCISSORS

**SofTouch Suture Scissors 17cm Castroviejo
Round Handle Curved Blunt Tips Duck Billed**

SKU: N-4CXF-7

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**SofTouch Suture Scissors 13cm Castroviejo
Flat Handle Curved Blunt Tips Duck Bill**

SKU: N-4CXF

Enquire through EndoTech NZ



MICROSURGICAL SCISSORS

**13cm Scissors Castroviejo Flat Handle
Curved Sharp Tips**

SKU: N-4C

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Micro forceps

1-6 OF 8 PRODUCTS



MICRO FORCEPS

15cm Micro Forceps ThumLok® Carbide Flat Handle 45 deg N/S

SKU: 45SL/M

Enquire through EndoTech NZ



MICRO FORCEPS

15cm Micro Forceps for Inlays/Onlays Carbide Flat Handle 45 deg N/S

SKU: 45-SM

Enquire through EndoTech NZ



MICRO FORCEPS

16cm Universal Forceps ThumLok® Carbide Flat Handle 45 deg E/W

SKU: 45SPL/M

Enquire through EndoTech NZ



MICRO FORCEPS

14.5cm Universal Corn Forceps Castroviejo Lock Diamond Flat-Handle

SKU: CORN

Enquire through EndoTech NZ



MICRO FORCEPS

17cm Universal Corn Forceps Castroviejo Lock Diamond Flat-Handle 45 deg Curved

SKU: CORN/45

Enquire through EndoTech NZ



MICRO FORCEPS

15cm Universal Corn Forceps Castroviejo Lock Diamond Flat-Handle 45 deg Angled


SKU: CORN/45B

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Micro forceps continued

7-8 OF 8 PRODUCTS




MICRO FORCEPS

**18cm Lasner Tissue Forceps ThumLok®
Round Handle Straight, One Pin**

SKU: MA-1S

Enquire through EndoTech NZ



MICRO FORCEPS

**17cm Lasner Tissue Forceps ThumLok®
Round Handle 45°N/S, Double Prong**

SKU: MA-4S

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Accessories

1-6 OF 7 PRODUCTS



ACCESSORIES

11.9cm Retractor Straight

SKU: GRH

Enquire through EndoTech NZ



ACCESSORIES

11.6cm Retractor Angled

SKU: GRH-A

Enquire through EndoTech NZ



ACCESSORIES

15.5cm Suture Tying Forceps ThumLok® Carbide Round Handle

SKU: PLAF/R/F-L

Enquire through EndoTech NZ



ACCESSORIES

15.5cm Suture Tying Forceps Carbide Round Handle

SKU: PLAF/R/F-6

Enquire through EndoTech NZ



ACCESSORIES

18cm Suture Tying Forceps Carbide Round Handle

SKU: PLAF/R/F

Enquire through EndoTech NZ



ACCESSORIES

Chromic Gut Suture 5-0, Reverse Cutting, ..

SKU: V-GUT-5-0 Chromic

Enquire through EndoTech NZ

FULL ENDODONTIC RANGE

Accessories continued

7-7 OF 7 PRODUCTS



ACCESSORIES

Chromic Gut Suture 4-0, Reverse Cutting, C-6 Needle/19mm, 27 / 70cm..Absorbable Surgical Suture

SKU: V-GUT-4-0 Chromic

Enquire through EndoTech NZ

TECHNIQUE VIDEOS

Clickable video links for product familiarisation.

These links are included for clinician reference and product orientation. Each link opens the relevant YouTube video.

Endodontic file control



Dr. Marc Habib - Laschal EF-1-7 tactile file forceps

Related SKU: EF-1-7
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=O14DwrteW9M>



Laschal endodontic file forceps

Related SKU: 90AHF/L
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=D0TvF82CMR4>

Separated file retrieval / FXP



FXP troughing system

Related SKU: FXP/S
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=afvIQxp7Ls>



Separated file removal with the Laschal FXP system

Related SKU: FXP/S
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=jEdIVL3Mzsk>



Dr Ortiz Hughes demonstrates FXP retrieval

Related SKU: FXP/S
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=249Xa1SXE8>



Steiglitz retrieval forceps support

Related SKU: D-45SL/M
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=iJ9I7UQPtVs>

Microsurgical closure / suturing



SofTouch suture scissors

Related SKU: N-4XF
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=tLkcUML19kM>




Snagless needle holder

Related SKU: 7-TCLR/TL
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=jcQq5E0I4wl>



ThumLok® needle holder

Related SKU: 7-TCLR/TL
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=7VOFwNafd3Y>



Cutting Edge needle holder

Related SKU: 6-10RL
Laschal Surgical, LLC.
<https://www.youtube.com/watch?v=GOS4E0fQ0ZM>

TECHNIQUE VIDEOS

Additional product video links.

These links are included for clinician reference and product orientation. Each link opens the relevant YouTube video.

Surgical access / perio and implant support



Periotomes

Laschal Surgical, LLC.

https://www.youtube.com/watch?v=_iE2rTe9KT0



CORN forceps

Related SKU: CORN/45B

Laschal Surgical, LLC.

<https://www.youtube.com/watch?v=4rawUzvBKjk>



Tissue graft forceps

Related SKU: MA-15

Laschal Surgical, LLC.

https://www.youtube.com/watch?v=_kZobKokHGQ



ENDOTECH NZ SUPPORT

Talk to EndoTech NZ about Laschal instruments.

Tell us the procedure you are preparing for and the access challenge you want to solve.
We can help identify suitable Laschal instruments and confirm availability.

[Contact EndoTech NZ](#)

[Download catalogue](#)